

See discussions, stats, and author profiles for this publication at: <https://www.researchgate.net/publication/341830560>

AutoDock Tutorial

Method · June 2020
DOI: 10.13140/RG.2.2.27339.21289

CITATIONS
0

READS
7,454

1 author:



[Ammar Elmezayen](#)
University of Texas Southwestern Medical Center
11 PUBLICATIONS 505 CITATIONS

SEE PROFILE

Some of the authors of this publication are also working on these related projects:



The discovery of novel class IIa HDACs inhibitors [View project](#)



Microsomal prostaglandin E synthase-1 (mPGES-1) or Prostaglandin E synthase Inhibitors [View project](#)

AUTODOCK TUTORIAL

Prepared by
Ammar Elmezayen

AUTODOCK TUTORIAL

- *At the start of this tutorial, all you need is a directory (e.g., “Docking”) which consists of two files; **protein.pdb** and **ligand.mol2**.*
- 1. **Open** ADT tool of AutoDock.
- 2. **File -> Preferences -> Set:** change the startup directory.
- 3. **Ligand -> Input -> Open:** open your ligand mol2 file.
- 4. **Ligand -> Torsion Tree -> Choose Torsions:** decide on the rotatable bonds to be considered.
- 5. **Ligand -> Output -> Save as PDBQT:** save your ligand file with “pdbqt” extension (e.g. ligand.pdbqt). Type in the extension manually.

AUTODOCK TUTORIAL

- **6. Grid -> Macromolecule -> Open:** Open your protein pdb file (e.g. protein.pdb), and ADT tool will prompt a new window to save your protein into pdbqt file.
- **7. Grid -> Set Map Types -> Choose Ligand:** select the ligand molecule from the list shown.
- **8. Grid -> Grid Box:** Select the number of grid points in x, y and z directions. Usually, the default 40 is taken. However, you can increase or decrease that number in any direction. From “**Center**” menu, select “**Center on ligand**”. Leave other parameters on default.

(Observe your grid box. Make sure that it covers the whole ligand as well as the binding site).

AUTODOCK TUTORIAL


- **File -> Close saving current:** save your settings before closing. Otherwise, your settings will be lost if the program is closed.
- **10. Grid -> Output -> Save GPF:** save your gpf file (e.g. docking.gpf). Type in the extension manually.
- **11. Check out** the content of Grid Parameter File, docking.gpf (by Notepad).
- **12. Docking -> Macromolecule -> Set Rigid Filename:** Select your macromolecule protein file (e.g. protein.pdbqt).
- **13. Docking -> Ligand -> Choose:** select your ligand. Leave everything on default and select **Accept**.

AUTODOCK TUTORIAL

- 14. **Docking -> Search Parameters -> Genetic Algorithm:** set the number of GA Runs as 20.
- Also, set the maximum number of “evals” to *long = 25000000* depending on the number rotatable bonds. Leave other parameters on default and select **Accept**.
- 15. **Docking -> Docking Parameters:** Leave on default.
- 16. **Docking -> Output -> Lamarckian GA:** save your dpf file (e.g. docking.dpf). Type in the extension manually.
- 17. **Close** your ADT tool.

AUTODOCK TUTORIAL

- In your working directory you should now have the following four files:
 - ligand.pdbqt • receptor.pdbqt • docking.gpf • docking.dpf
- 18. **Open “Windows PowerShell”:** **File -> Open Windows PowerShell -> Open Windows PowerShell as administrator.**



AUTODOCK TUTORIAL

- While you are inside your working directory, type the following Win/unix command:

Windows:

`.\autogrid4.exe -p docking.gpf -l docking.glg`

Unix:

`autogrid4 -p docking.gpf -l docking.glg`

- Once the program is completed, new files should appear in your working directory including the map files.

AUTODOCK TUTORIAL

- These grid map files are necessary for the docking run. Now you are ready to dock. Type the following:

Windows:

```
.\autodock4.exe -p docking.dpf -l docking.dlg
```

Unix:

```
autodock4 -p docking.dpf -l docking.dlg
```

HAPPY DOCKING

